

Nur für Forschungszwecke

# TIRAP Monoklonaler Antikörper

Katalog-Nr.:60016-1-Ig



## Allgemeine Informationen

<b>Katalog-Nr.:</b> 60016-1-Ig	<b>GenBank-Zugangsnummer:</b> BC032474	<b>Reinigungsmethode:</b> Protein-G-Reinigung
<b>Größe:</b> 150ul , Konzentration: 1000 µg/ml durch die Bradford-Methode mit BSA als Standard;	<b>GeneID (NCBI):</b> 114609	<b>CloneNo.:</b> 4H5E12
<b>Wirt:</b> Maus	<b>Vollständiger Name:</b> toll-interleukin 1 receptor (TIR) domain containing adaptor protein	
<b>Isotyp:</b> IgG1	<b>Berechnete Masse:</b> 256 aa, 28 kDa	
<b>Immunogen Katalognummer:</b> AG0771	<b>Beobachtete Masse:</b> 35 kDa	

## Anwendungen

**Geprüfte Anwendungen:**  
ELISA

**Getestete Reaktivität:**  
Human

## Hintergrundinformationen

### Lagerung

**Lagerungsbedingungen:**  
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

**Lagerungspuffer:**  
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

