

Nur für Forschungszwecke

RARRES2 Monoklonaler Antikörper

Katalog-Nr.: 60018-1-Ig

1 Publikationen



Allgemeine Informationen

Katalog-Nr.:

60018-1-Ig

Größe:

150ul

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG0261

GenBank-Zugangsnummer:

BC000069

GeneID (NCBI):

5919

Vollständiger Name:

retinoic acid receptor responder
(tazarotene induced) 2

Berechnete Masse:

163 aa, 19 kDa

Reinigungsmethode:

Caprylsäure/Ammoniumsulfat-
Präzipitation

CloneNo.:

2C4D9

Anwendungen

Geprüfte Anwendungen:

ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human

Hintergrundinformationen

Bemerkenswerte Veröffentlichungen

Verfasser

Xiao-Hui Song

Pubmed ID

34957549

Journal

J Cell Physiol

Anwendung

WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

