

Nur für Forschungszwecke

PRIM1 Monoklonaler Antikörper

Katalog-Nr.: 60021-1-Ig



Allgemeine Informationen

Katalog-Nr.:

60021-1-Ig

Größe:

150ul

Wirt:

Maus

Isotyp:

IgG2b

Immunogen Katalognummer:

AG1124

GenBank-Zugangsnummer:

BC005266

GeneID (NCBI):

5557

Vollständiger Name:

primase, DNA, polypeptide 1 (49kDa)

Berechnete Masse:

50 kDa

Reinigungsmethode:

Caprylsäure/Ammoniumsulfat-
Präzipitation

CloneNo.:

5H2C8

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

