

Allgemeine Informationen

Katalog-Nr.: 60033-1-Ig	GenBank-Zugangsnummer: BE284390	Reinigungsmethode: Protein-G-Reinigung
Größe: 150ul, Konzentration: 1000 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 289094	CloneNo.: 4D3E10
Wirt: Maus	Vollständiger Name: hemicentin 1	
Isotyp: IgG1	Berechnete Masse: 612 kDa	
Immunogen Katalognummer: AG1770	Beobachtete Masse: 58-60 kDa	

Anwendungen

Geprüfte Anwendungen:
ELISA

In Publikationen genannte Anwendungen:
FC, IHC

Getestete Reaktivität:
Human, Maus

Zitierte Arten:
Human

Hintergrundinformationen

This gene encodes a large extracellular member of the immunoglobulin superfamily. This protein were found in ECM of many epithelial cells, including embryonic trophoctoderm and adult skin and tongue. It may function in in many processes such as stabilization of the germline syncytium, anchorage of mechanosensory neurons to the epidermis, and organization of hemidesmosomes in the epidermis.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Margherita Sisto	19190085	Int Immunol	FC
Kortvely Elod E	19696174	Invest Ophthalmol Vis Sci	IHC

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

