

Nur für Forschungszwecke

Hmcn1 Monoklonaler Antikörper

Katalog-Nr.: 60033-2-Ig



Allgemeine Informationen

Katalog-Nr.:

60033-2-Ig

Größe:

150ul

Wirt:

Maus

Isotyp:

IgG2a

Immunogen Katalognummer:

AG1770

GenBank-Zugangsnummer:

BE284390

GeneID (NCBI):

289094

Vollständiger Name:

hemicentin 1

Berechnete Masse:

612 kDa

Reinigungsmethode:

Caprylsäure/Ammoniumsulfat-
Präzipitation

CloneNo.:

5F9F4

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

This gene encodes a large extracellular member of the immunoglobulin superfamily. This protein were found in ECM of many epithelial cells, including embryonic trophectoderm and adult skin and tongue. It may function in many processes such as stabilization of the germline syncytium, anchorage of mechanosensory neurons to the epidermis, and organization of hemidesmosomes in the epidermis.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

