

Nur für Forschungszwecke

# Hmcn1 Monoklonaler Antikörper

Katalog-Nr.:60033-3-Ig



## Allgemeine Informationen

Katalog-Nr.:  
60033-3-Ig

Größe:  
150ul, Konzentration: 1000 µg/ml  
durch die Bradford-Methode mit BSA  
als Standard;

Wirt:  
Maus

Isotyp:  
IgM

Immunogen Katalognummer:  
AG1770

GenBank-Zugangsnummer:  
BE284390

GeneID (NCBI):  
289094

Vollständiger Name:  
hemicentin 1

Berechnete Masse:  
612 kDa

Reinigungsmethode:

Caprylsäure/Ammoniumsulfat-  
Präzipitation

CloneNo.:  
2E8E4

## Anwendungen

Geprüfte Anwendungen:  
ELISA

Getestete Reaktivität:  
Human

## Hintergrundinformationen

This gene encodes a large extracellular member of the immunoglobulin superfamily. This protein were found in ECM of many epithelial cells, including embryonic trophectoderm and adult skin and tongue. It may function in in many processes such as stabilization of the germline syncytium, anchorage of mechanosensory neurons to the epidermis, and organization of hemidesmosomes in the epidermis.

## Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

