

Nur für Forschungszwecke

MAGEC2 Monoklonaler Antikörper

Katalog-Nr.:60040-1-Ig



Allgemeine Informationen

Katalog-Nr.:

60040-1-Ig

Größe:

150ul , Konzentration: 1000 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Maus

Isotyp:

IgG2a

Immunogen Katalognummer:

AG0418

GenBank-Zugangsnummer:

BC005891

GeneID (NCBI):

51438

Vollständiger Name:

melanoma antigen family C, 2

Berechnete Masse:

373 aa, 41 kDa

Reinigungsmethode:

Caprylsäure/Ammoniumsulfat-Präzipitation

CloneNo.:

2H5A12

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

Melanoma antigen family C, 2(MAGEC2, synonyms: CT10, HCA587, MAGEE1, MGC13377) is related to members of the MAGEC gene family. MAGEC2 is not expressed in normal tissues, except for testis, and is expressed in tumors of various histological types. This gene and other MAGEC genes are clustered on chromosome Xq26-q27.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

