

Nur für Forschungszwecke

PSMG1 Monoklonaler Antikörper

Katalog-Nr.: 60048-1-Ig



Allgemeine Informationen

Katalog-Nr.:

60048-1-Ig

Größe:

150ul, Konzentration: 1000 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Maus

Isotyp:

IgM

Immunogen Katalognummer:

AG0303

GenBank-Zugangsnummer:

BC003619

GeneID (NCBI):

8624

Vollständiger Name:

proteasome (prosome, macropain)
assembly chaperone 1

Berechnete Masse:

288 aa, 33 kDa

Beobachtete Masse:

44 kDa

Reinigungsmethode:

Caprylsäure/Ammoniumsulfat-
Präzipitation

CloneNo.:

1B6E1

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

