

Nur für Forschungszwecke

Rab18 Monoklonaler Antikörper

Katalog-Nr.: 60057-1-Ig **4 Publikationen**



Allgemeine Informationen

Katalog-Nr.: 60057-1-Ig	GenBank-Zugangsnummer: BC015014	Reinigungsmethode: Caprylsäure/Ammoniumsulfat-Präzipitation
Größe: 150ul, Konzentration: 1000 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 22931	CloneNo.: 4D5B8
Wirt: Maus	Vollständiger Name: RAB18, member RAS oncogene family	
Isotyp: IgG1	Berechnete Masse: 206 aa, 23 kDa	
Immunogen Katalognummer: AG1841		

Anwendungen

Geprüfte Anwendungen:
ELISA

In Publikationen genannte Anwendungen:
IF, WB

Getestete Reaktivität:
Human, Maus

Zitierte Arten:
Affe, Maus

Hintergrundinformationen

RAB18 is a member of a family of Ras-related small GTPases. RAB18 plays a role in apical endocytosis/recycling, and may be implicated in transport between the plasma membrane and early endosomes. It also plays a key role in eye and brain development and neurodegeneration. Defects in RAB18 are associated with Warburg micro syndrome type 3.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Yu Yan	33045621	Redox Biol	IF
Andreas Gerondopoulos	24891604	J Cell Biol	WB
Dijin Xu	29367353	J Cell Biol	IF

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

