

Nur für Forschungszwecke

LAMR1,RPSA Monoklonaler Antikörper

Katalog-Nr.:60068-1-Ig



Allgemeine Informationen

Katalog-Nr.:
60068-1-Ig

Größe:
150ul , Konzentration: 807 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Maus

Isotyp:
IgG2b

Immunogen Katalognummer:
AG1055

GenBank-Zugangsnummer:
BC008867

GeneID (NCBI):
3921

Vollständiger Name:
ribosomal protein SA

Berechnete Masse:
33 kDa

Beobachtete Masse:
40 kDa

Reinigungsmethode:

Caprylsäure/Ammoniumsulfat-
Präzipitation

CloneNo.:
5H3G8

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human, Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

