

Nur für Forschungszwecke

PFN2 Monoklonaler Antikörper

Katalog-Nr.:60094-1-Ig **3 Publikationen**



Allgemeine Informationen

Katalog-Nr.: 60094-1-Ig	GenBank-Zugangsnummer: BC018049	Reinigungsmethode: Protein-A-Reinigung
Größe: 150ul	GeneID (NCBI): 5217	CloneNo.: 1D9C4
Wirt: Maus	Vollständiger Name: profilin 2	
Isotyp: IgM	Berechnete Masse: 140 aa, 15 kDa	
Immunogen Katalognummer: AG5043		

Anwendungen

Geprüfte Anwendungen:
ELISA

In Publikationen genannte Anwendungen:
IF, WB

Getestete Reaktivität:
Human, Maus

Zitierte Arten:
Human, Maus

Hintergrundinformationen

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Min Yu	29773903	Oncogene	WB
Sonia Afroz	27136678	Elife	IF
Xiaobin Cui	29238898	Mol Cell Biochem	WB

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

