

Nur für Forschungszwecke

APOBEC3G Monoklonaler Antikörper

Katalog-Nr.:60100-2-Ig



Allgemeine Informationen

Katalog-Nr.: 60100-2-Ig	GenBank-Zugangsnummer: BC024268	Reinigungsmethode: Protein-A-Reinigung
Größe: 150ul	GeneID (NCBI): 60489	CloneNo.: 5G10A8
Wirt: Maus	Vollständiger Name: apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3G	
Isotyp: IgG2b	Berechnete Masse: 46 kDa	
Immunogen Katalognummer: AG7770		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

APOBEC3G (apolipoprotein B mRNA-editing enzyme, catalytic polypeptide-like 3G, also known as CEM15), a member of the APOBEC family of cytidine deaminase-editing enzymes, is a potent DNA mutator that can be packaged into budding retroviruses. HA-tagged APOBEC3G can be expressed in 293T cells, which do not express APOBEC3G endogenously, in addition 22- and 24-kDa proteolytic degradation fragments can be detected by western blot (PMID:14672928). It has 2 isoforms produced by alternative splicing.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

