

Nur für Forschungszwecke

TWEAKR Monoklonaler Antikörper

Katalog-Nr.:60104-1-Ig



Allgemeine Informationen

Katalog-Nr.: 60104-1-Ig	GenBank-Zugangsnummer: BC002718	Reinigungsmethode: Caprylsäure/Ammoniumsulfat- Präzipitation
Größe: 150ul	GeneID (NCBI): 51330	CloneNo.: 2E10C10
Wirt: Maus	Vollständiger Name: tumor necrosis factor receptor superfamily, member 12A	
Isotyp: IgG2b	Berechnete Masse: 129 aa, 14 kDa	
Immunogen Katalognummer: AG2689		

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



human heart tissue were subjected to SDS PAGE followed by western blot with 60104-1-Ig (TWEAKR antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.