

Nur für Forschungszwecke

# FXVD6 Monoklonaler Antikörper

Katalog-Nr.:60111-1-Ig



## Allgemeine Informationen

<b>Katalog-Nr.:</b> 60111-1-Ig	<b>GenBank-Zugangsnummer:</b> BC018652	<b>Reinigungsmethode:</b> Caprylsäure/Ammoniumsulfat- Präzipitation
<b>Größe:</b> 150ul	<b>GeneID (NCBI):</b> 53826	<b>CloneNo.:</b> 6H10F4
<b>Wirt:</b> Maus	<b>Vollständiger Name:</b> FXVD domain containing ion transport regulator 6	
<b>Isotyp:</b> IgG2a	<b>Berechnete Masse:</b> 95 aa, 11 kDa	
<b>Immunogen Katalognummer:</b> AG8538		

## Anwendungen

**Geprüfte Anwendungen:**  
ELISA  
**Getestete Reaktivität:**  
Human

## Hintergrundinformationen

### Lagerung

**Lagerungsbedingungen:**  
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil  
**Lagerungspuffer:**  
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.  
Aliquotieren ist nicht notwendig bei -20°C Lagerung

**\*\*\* 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com  
W: ptglab.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

