

Nur für Forschungszwecke

Hepsin Monoklonaler Antikörper

Katalog-Nr.: 60114-1-Ig



Allgemeine Informationen

Katalog-Nr.:

60114-1-Ig

Größe:

150ul, Konzentration: 1000 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Maus

Isotyp:

IgG1

GenBank-Zugangsnummer:

NM_002151

GeneID (NCBI):

3249

Vollständiger Name:

hepsin

Berechnete Masse:

45 kDa

Beobachtete Masse:

51 kDa

Reinigungsmethode:

Caprylsäure/Ammoniumsulfat-
Präzipitation

CloneNo.:

2D3A9

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

Hepsin, a type II membrane-associated protein, is a cell surface proteases. It acts directly by degrading the extracellular matrix of surrounding cells and also functions as an activator of other extracellular proteases. It was reported to be involved in multiple biological processes.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

