

Nur für Forschungszwecke

Kv1.2 Monoklonaler Antikörper

Katalog-Nr.:60120-1-Ig



Allgemeine Informationen

Katalog-Nr.:
60120-1-Ig

Größe:
150ul, Konzentration: 500 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:
Maus

Isotyp:
IgG1

Immunogen Katalognummer:
AG5477

GenBank-Zugangsnummer:
BC043564

GeneID (NCBI):
3737

Vollständiger Name:
potassium voltage-gated channel,
shaker-related subfamily, member 2

Berechnete Masse:
499 aa, 57 kDa

Beobachtete Masse:
75 kDa

Reinigungsmethode:
Caprylsäure/Ammoniumsulfat-
Präzipitation

CloneNo.:
1D2D11

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

