

Nur für Forschungszwecke

IL13RA2,CD213A2 Monoklonaler Antikörper

Katalog-Nr.:60140-1-Ig



Allgemeine Informationen

Katalog-Nr.:

60140-1-Ig

Größe:

150ul , Konzentration: 500 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG1534

GenBank-Zugangsnummer:

BC020739

GeneID (NCBI):

3598

Vollständiger Name:

interleukin 13 receptor, alpha 2

Berechnete Masse:

380 aa, 44 kDa

Reinigungsmethode:

Caprylsäure/Ammoniumsulfat-
Präzipitation

CloneNo.:

2F10B3

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

IL13RA2, also named as IL13R and CD213a2, belongs to the type I cytokine receptor family and Type 5 subfamily. It binds as a monomer with high affinity to interleukin-13 (IL13), but not to interleukin-4 (IL4). The IL13RA2 gene is often overexpressed in brain tumors (PMID:19895199).

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

