

Nur für Forschungszwecke

KCNH7 Monoklonaler Antikörper

Katalog-Nr.: 60166-1-Ig



Allgemeine Informationen

Katalog-Nr.:

60166-1-Ig

Größe:

150ul, Konzentration: 500 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG4723

GenBank-Zugangsnummer:

BC035815

GeneID (NCBI):

90134

Vollständiger Name:

potassium voltage-gated channel,
subfamily H (eag-related), member 7

Berechnete Masse:

1196 aa, 135 kDa

Beobachtete Masse:

120-130 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

2F2G1

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

