

Nur für Forschungszwecke

GCDFP-15/PIP Monoklonaler Antikörper



Katalog-Nr.: 60210-1-Ig

1 Publikationen

Allgemeine Informationen

Katalog-Nr.:

60210-1-Ig

Größe:

150ul, Konzentration: 547 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG8980

GenBank-Zugangsnummer:

BC010950

GeneID (NCBI):

5304

Vollständiger Name:

prolactin-induced protein

Berechnete Masse:

146 aa, 17 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

4F5C11

Anwendungen

Geprüfte Anwendungen:

ELISA

In Publikationen genannte Anwendungen:

IF

Getestete Reaktivität:

Human, Maus

Zitierte Arten:

Rind

Hintergrundinformationen

GCDFP-15 (gross cystic disease fluid protein 15), also known as PIP (prolactin-induced protein) is a secretory glycoprotein expressed in benign and malignant breast tumor tissues and in some normal exocrine organs such as sweat, salivary, and lacrimal glands (PMID: 12393800). GCDFP-15 expression is increased by prolactin and androgen and inhibited by estrogen (PMID: 18854942). The expression is also regulated by interleukins. GCDFP-15 is a marker of apocrine differentiation and is frequently used for assessment of metastases or regional recurrences of breast origin.

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
S L Pratt	25533929	Theriogenology	IF

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

