

Nur für Forschungszwecke

CPS1 Monoklonaler Antikörper

Katalog-Nr.:65011-1-Ig **1 Publikationen**



Allgemeine Informationen

Katalog-Nr.: 65011-1-Ig	GenBank-Zugangsnummer: BC020695	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , 667 µg/ml	GeneID (NCBI): 1373	Empfohlene Verdünnungen: WB 1:200-1:1000
Wirt: Maus	Vollständiger Name: carbamoyl-phosphate synthetase 1, mitochondrial	
Isotyp: IgG	Berechnete Masse: 165 kDa	

Anwendungen

Geprüfte Anwendungen: WB,ELISA	Positivkontrollen: WB : Human liver,
In Publikationen genannte Anwendungen: WB	
Getestete Reaktivität: Human	
Zitierte Arten: Human	

Hintergrundinformationen

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Jana Deitersen	34059630	Cell Death Dis	WB

Lagerung

Lagerungsbedingungen:
Bei 2-8°C lagern.
Lagerungspuffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

