

Nur für Forschungszwecke

APOA5 Monoklonaler Antikörper

Katalog-Nr.:65015-1-Ig



Allgemeine Informationen

Katalog-Nr.: 65015-1-Ig	GenBank-Zugangsnummer: BC101787	Reinigungsmethode: Antigen-Affinitätsreinigung
Größe: 150ul , 500 µg/ml	GeneID (NCBI): 116519	Empfohlene Verdünnungen: WB 1:500-1:2000
Wirt: Maus	Vollständiger Name: apolipoprotein A-V	
Isotyp: IgG1	Berechnete Masse: 366 aa, 41 kDa	
	Beobachtete Masse: 39 kDa-41 kDa	

Anwendungen

Geprüfte Anwendungen: WB,ELISA	Positivkontrollen: WB : Human serum,
Getestete Reaktivität: Human	

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei 2-8°C lagern.
Lagerungspuffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

