

DR5 Monoklonaler Antikörper

Katalog-Nr.: 65046-1-Ig

Allgemeine Informationen

Katalog-Nr.:

65046-1-Ig

Größe:

150 µl, 1000 µg/ml

Wirt:

Maus

Isotyp:

IgG1

GenBank-Zugangsnummer:

BC001281

GeneID (NCBI):

8795

Vollständiger Name:

tumor necrosis factor receptor
superfamily, member 10b

Berechnete Masse:

48 kDa

Beobachtete Masse:

58 kDa

Reinigungsmethode:

Protein-G-Reinigung

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

TNFRSF10B, also known as CD262, DR5, TRAILR2, TRICK2 and KILLER, is a widely expressed single-pass type I membrane protein belonging to the tumour necrosis factor receptor superfamily (TNFRSF). It is a receptor for TNF-related apoptosis-inducing ligand (TRAIL), which is a member of the tumor necrosis factor (TNF) family of cytokines and induces apoptosis in a wide variety of cells (PMID: 9311998). TNFRSF10B contains two extracellular cysteine-rich repeats, typical for TNF receptor (TNFR) family members, and a cytoplasmic death domain (DD), through which TNFRSF10B is capable to transmit the apoptotic signal (PMID: 9311998; 20531300).

Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern.

Lagerungspuffer:

0.1M Tris-HCl, 0.15M NaCl, 0.04% merthiolate and 50% glycerol pH 7.4.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

