

Nur für Forschungszwecke

CD107a / LAMP1 Monoklonaler Antikörper



Katalog-Nr.: 65050-1-Ig **1 Publikationen**

Allgemeine Informationen

Katalog-Nr.: 65050-1-Ig	GenBank-Zugangsnummer: BC006785	Reinigungsmethode: Affinitätsreinigung
Größe: 500ug, 0.5 mg/ml	GeneID (NCBI): 16783	CloneNo.: 1D4B
Wirt: Ratte	Vollständiger Name: lysosomal-associated membrane protein 1	Empfohlene Verdünnungen: IF 1:50-1:500
Isotyp: IgG2a		

Anwendungen

Geprüfte Anwendungen: FC, IF	Positivkontrollen: IF : NIH/3T3-Zellen,
In Publikationen genannte Anwendungen: IF	
Getestete Reaktivität: Maus	
Zitierte Arten: Maus	

Hintergrundinformationen

LAMP1 (CD107a) is a heavily glycosylated membrane protein enriched in the lysosomal membrane. LAMP1 is extensively glycosylated with asparagine-linked oligosaccharides which protect it from intracellular proteolysis (PMID: 10521503). Although LAMP1 is expressed largely in the endosome-lysosomal membrane of cells, it is also found on the plasma membrane (PMID: 16168398). Elevated LAMP1 expression at the cell surface has also been detected during platelet and granulocytic cell activation, as well as in some tumor cells (PMID: 29085473). LAMP1 functions to provide selectins with carbohydrate ligands. This protein has also been shown to be a marker of degranulation on lymphocytes such as CD8+ and NK cells and may also play a role in tumor cell differentiation and metastasis (PMID: 18835598; 29085473; 9426697).

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Rui Hua	36924084	J Cell Physiol	IF



Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Stable for one year after shipment.
Lagerungspuffer:
Borat-gepufferte Kochsalzlösung, pH 8,2.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ NIH/3T3 cells were intracellularly stained with 0.50 ug Anti-Mouse CD107a (65050-1-Ig, Clone: 1D4B) (red) or 0.50 ug isotype control antibody (blue), and Fluorescein (FITC)-conjugated Affinipure Goat Anti-Rat IgG(H+L) with dilution of 1:1000. Cells were fixed with 90% MeOH.

Immunofluorescent analysis of (4% PFA) fixed NIH/3T3 cells using 65050-1-Ig (CD107a antibody) at dilution of 1:100 and Rhodamine (TRITC)-conjugated Goat Anti-Rat IgG(H+L).