

Nur für Forschungszwecke

CD115 (c-fms) Monoklonaler Antikörper

Katalog-Nr.:65095-1-Ig



Allgemeine Informationen

Katalog-Nr.:

65095-1-Ig

Größe:

100ug , 500 µg/ml

Wirt:

Ratte

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC036343

GeneID (NCBI):

12978

Vollständiger Name:

colony stimulating factor 1 receptor

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

AFS98

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Stable for one year after shipment.

Lagerungspuffer:

PBS mit 0,1% Natriumazid, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1×10^6 Balb/c mouse peritoneal macrophages were surface stained with 0.5 μ g Anti-Mouse CD115 (c-fms) (65095-1-Ig, Clone: AFS98) and FITC anti-rat IgG2a Antibody at dilution 1:200 (green), or stained with 0.5 μ g isotype control and FITC anti-rat IgG2a Antibody at dilution 1:200 (black). Cells were not fixed.