

Nur für Forschungszwecke

# CD19 Monoklonaler Antikörper

Katalog-Nr.: 65110-1-Ig



## Allgemeine Informationen

<b>Katalog-Nr.:</b> 65110-1-Ig	<b>GenBank-Zugangsnummer:</b> BC006338	<b>Reinigungsmethode:</b> Affinitätsreinigung
<b>Größe:</b> 100µg, 500 µg/ml	<b>GeneID (NCBI):</b> 930	<b>CloneNo.:</b> HIB19
<b>Wirt:</b> Maus	<b>Vollständiger Name:</b> CD19 molecule	
<b>Isotyp:</b> Mouse IgG1, kappa	<b>Berechnete Masse:</b> 556 aa, 61 kDa	

## Anwendungen

**Geprüfte Anwendungen:**  
FC

**Getestete Reaktivität:**  
Human

## Hintergrundinformationen

CD19 is a 95 kDa type I transmembrane glycoprotein belonging to the immunoglobulin superfamily (PMID: 2472450). It is expressed by B cells and follicular dendritic cells. CD19 is up-regulated at the step of B-lineage commitment during the differentiation of the hematopoietic stem cell, it remains on during subsequent stages of differentiation until finally down-regulated during terminal differentiation into plasma cells (PMID: 8528044). CD19 is involved in B cell development, activation and differentiation. It is the dominant component for the signaling complex on B cells that includes CD21 (CR2), CD81 (TAPA-1) and CD225 and acts as a critical co-receptor for BCR signal transduction (PMID: 23210908).

## Lagerung

**Lagerungsbedingungen:**  
Bei 2-8°C lagern.

**Lagerungspuffer:**  
PBS mit 0,1% Natriumazid, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> human peripheral blood lymphocytes were surface stained with 0.5 ug Anti-Human CD19 (65110-1-Ig, Clone: H1B19) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000 (green), or stained with 0.5 ug isotype control and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000 (black). Cells were not fixed.