

Nur für Forschungszwecke

CD24 Monoklonaler Antikörper

Katalog-Nr.:65121-1-Ig



Allgemeine Informationen

Katalog-Nr.: 65121-1-Ig	GenBank-Zugangsnummer: BC075622	Reinigungsmethode: Affinitätsreinigung
Größe: 100ug , 500 µg/ml	GeneID (NCBI): 12484	CloneNo.: M1/69
Wirt: Ratte	Vollständiger Name: CD24a antigen	
Isotyp: IgG2b, kappa		

Anwendungen

Geprüfte Anwendungen:
FC

Getestete Reaktivität:
Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Store at 2-8°C. Stable for one year after shipment.

Lagerungspuffer:
PBS mit 0,1% Natriumazid, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ BALB/c mouse splenocytes were surface stained with 0.5 ug Anti-Mouse CD24 (65121-1-Ig, Clone: M1/69) and FITC anti-rat IgG2b Antibody at dilution 1:200 (green), or stained with 0.5 ug isotype control and FITC anti-rat IgG2b Antibody at dilution 1:200 (black). Cells were not fixed.