

Nur für Forschungszwecke

CD116 Monoklonaler Antikörper

Katalog-Nr.:65160-1-Ig



Allgemeine Informationen

Katalog-Nr.:

65160-1-Ig

Größe:

100ug , 1 mg/ml

Wirt:

Maus

Isotyp:

Mouse IgG1 kappa

GenBank-Zugangsnummer:

BC002635

GeneID (NCBI):

1438

Vollständiger Name:

colony stimulating factor 2 receptor,
alpha, low-affinity (granulocyte-
macrophage)

Berechnete Masse:

46 kDa

Reinigungsmethode:

Mit Protein-A-

Affinitätschromatographie gereinigt

CloneNo.:

4H1

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern.

Lagerungspuffer:

TBS mit 0,1% Natriumazid.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human peripheral blood monocytes were surface stained with 0.5 ug Anti-Human CD116 (65160-1-Ig, Clone: 4H1) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000 (green), or stained with 0.5 ug isotype control and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000 (black). Cells were not fixed.