

Nur für Forschungszwecke

# CD1a Monoklonaler Antikörper

Katalog-Nr.:65172-1-Ig



## Allgemeine Informationen

Katalog-Nr.:

65172-1-Ig

Größe:

100ug, 1 mg/ml

Wirt:

Maus

Isotyp:

Mouse IgG1

GenBank-Zugangsnummer:

BC031645

GeneID (NCBI):

909

Vollständiger Name:

CD1a molecule

Berechnete Masse:

327 aa, 37 kDa

Reinigungsmethode:

Mit Protein-A-

Affinitätschromatographie gereinigt

CloneNo.:

HI149

## Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern.

Lagerungspuffer:

PBS mit 0,1% Natriumazid.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> MOLT-4 were surface stained with 0.5 ug Anti-Human CD1a (65172-1-Ig, Clone:HL149) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000 (green), or stained with 0.5 ug Mouse IgG1 isotype control and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000 (black). Cells were not fixed.