

Nur für Forschungszwecke

CD1c Monoklonaler Antikörper

Katalog-Nr.:65176-1-Ig



Allgemeine Informationen

Katalog-Nr.:

65176-1-Ig

Größe:

100ug, 1 mg/ml

Wirt:

Maus

Isotyp:

Mouse IgG1 kappa

GenBank-Zugangsnummer:

BC126465

GeneID (NCBI):

911

Vollständiger Name:

CD1c molecule

Reinigungsmethode:

Mit Protein-A-

Affinitätschromatographie gereinigt

CloneNo.:

L161

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei 2-8°C lagern.

Lagerungspuffer:

PBS mit 0,1% Natriumazid.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human peripheral blood lymphocytes were surface stained with APC Anti-Human CD19 (APC-65110, Clone: HIB19), and 1.5 ug Anti-Human CD1c (65176-1-Ig, Clone: L161) and FITC Anti-Mouse IgG1 Antibody. Cells were not fixed.



1X10⁶ human peripheral blood lymphocytes were surface stained with APC Anti-Human CD19 (APC-65110, Clone: HIB19), and 1.5 ug Mouse IgG1 Isotype Control and FITC Anti-Mouse IgG1 Antibody. Cells were not fixed.