

Nur für Forschungszwecke

CD41 Monoklonaler Antikörper

Katalog-Nr.:65189-1-Ig



Allgemeine Informationen

Katalog-Nr.:

65189-1-Ig

Größe:

100ug, 0.5 mg/ml

Wirt:

Ratte

Isotyp:

IgG1, kappa

GenBank-Zugangsnummer:

BC120493

GeneID (NCBI):

16399

Vollständiger Name:

integrin alpha 2b

Reinigungsmethode:

Mit Protein-G-Affinitätschromatographie aufgereinigt

CloneNo.:

MWReg30

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Maus

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ mouse splenocytes were surface stained with 0.5 μ g Anti-Mouse CD41 (65189-1-Ig, Clone:MWReg30) (red), or 0.5 μ g Isotype Control. Cells were not fixed.