

Nur für Forschungszwecke

IgG2a Isotype Control Monoklonaler Antikörper



Katalog-Nr.: 65208-1-Ig

Allgemeine Informationen

Katalog-Nr.: 65208-1-Ig	GenBank-Zugangsnummer: GeneID (NCBI):	Reinigungsmethode: Affinitätsreinigung
Größe: 500ug, 0.5 mg/mL	Vollständiger Name:	CloneNo.: C1.18.4
Wirt: Maus		
Isotyp: IgG2a, κ		

Anwendungen

Geprüfte Anwendungen:
FC
Getestete Reaktivität:

Hintergrundinformationen

The C1.18.4 immunoglobulin is useful as an isotype-matched control. The C.1.18.4 immunoglobulin has an unknown binding specificity and is used as an isotype control for mouse IgG2a antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Lagerung

Lagerungsbedingungen:
Bei 2-8°C lagern.
Lagerungspuffer:
Phosphatbasierter Puffer mit 0,09% Natriumazid, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human peripheral blood lymphocytes were surface stained with 0.5 ug Mouse IgG2a Isotype Control (65208-1-Ig, Clone: C1.18.4) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000 (black), or stained with 0.5 ug Anti-Human CD45 (65064-1-Ig, Clone: F10-89-4) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000 (green). Cells were not fixed.