

Nur für Forschungszwecke

# CD71 Monoklonaler Antikörper

Katalog-Nr.:65235-1-Ig



## Allgemeine Informationen

Katalog-Nr.:

65235-1-Ig

Größe:

100ug, 0.5 mg/ml

Wirt:

Maus

Isotyp:

IgG1, kappa

GenBank-Zugangsnummer:

BC001188

GeneID (NCBI):

7037

Vollständiger Name:

transferrin receptor (p90, CD71)

Berechnete Masse:

85 kDa

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

OKT9

## Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

## Hintergrundinformationen

### Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Ausgewählte Validierungsdaten



1X10<sup>6</sup> day 3 PHA treated human PBMCs were surface stained with 0.2 ug Anti-Human CD71 (65235-1-Ig, Clone:OKT9) and CL594-Goat anti-Mouse IgG at dilution 1:1000. Cells were not fixed. Lymphocytes were gated. Blue line indicates unstained stimulated cells. Black dashed line indicates unstimulated stained cells.