

Nur für Forschungszwecke

CD155 Monoklonaler Antikörper

Katalog-Nr.:65251-1-Ig



Allgemeine Informationen

Katalog-Nr.:

65251-1-Ig

Größe:

100ug, 0.5 mg/mL

Wirt:

Maus

Isotyp:

IgG1, kappa

GenBank-Zugangsnummer:

BC015542

GeneID (NCBI):

5817

Vollständiger Name:

poliovirus receptor

Berechnete Masse:

45 kDa

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

SKII.4

Anwendungen

Geprüfte Anwendungen:

FC, ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ U937 cells were surface stained with 0.2 ug Anti-Human CD155 (65251-1-Ig, Clone:SKII.4) and CoraLite®594-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000. Cells were not fixed.