

Nur für Forschungszwecke

CD26 Monoklonaler Antikörper

Katalog-Nr.:65267-1-Ig



Allgemeine Informationen

Katalog-Nr.:

65267-1-Ig

Größe:

100ug, 0.5 mg/mL

Wirt:

Maus

Isotyp:

IgG2a, kappa

GenBank-Zugangsnummer:

BC013329

GeneID (NCBI):

1803

Vollständiger Name:

dipeptidyl-peptidase 4

Berechnete Masse:

88 kDa

Reinigungsmethode:

Affinitätsreinigung

CloneNo.:

BA5b

Anwendungen

Geprüfte Anwendungen:

FC

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Store at 2-8°C. Stable for one year after shipment.

Lagerungspuffer:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



1X10⁶ human PBMCs were surface stained with 0.2 ug Anti-Human CD26 (65267-1-Ig, Clone:BA5b) and CoraLite®594-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000. Cells were not fixed. Lymphocytes were gated.



1X10⁶ human PBMCs were surface stained with 0.2 ug Anti-Human CD26 (65267-1-Ig, Clone:BA5b) and CoraLite®594-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000. Cells were not fixed. Lymphocytes were gated.