

Nur für Forschungszwecke

PAX4 Monoklonaler Antikörper

Katalog-Nr.:66064-2-Ig



Allgemeine Informationen

Katalog-Nr.: 66064-2-Ig	GenBank-Zugangsnummer: BC074761	Reinigungsmethode: Protein-A-Reinigung
Größe: 150ul , Konzentration: 873 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 5078	CloneNo.: 1D5C7
Wirt: Maus	Vollständiger Name: paired box 4	
Isotyp: IgG2b	Berechnete Masse: 343 aa, 37 kDa	
Immunogen Katalognummer: AG17007		

Anwendungen

Geprüfte Anwendungen:
ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

PAX4, a paired-box transcription factor, is a key regulator of pancreatic islet cell growth and differentiation. Differentiation of endoderm-derived endocrine pancreas is mediated through PAX4 and PAX6. PAX4 may act as a PAX6 repressor due to the competition for binding sites and lower transactivation potential of Pax-4. PAX4 has the potential to function as a tumor suppressor in human melanoma.

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

