

Nur für Forschungszwecke

IL-1F6 Monoklonaler Antikörper

Katalog-Nr.:66147-1-Ig



Allgemeine Informationen

Katalog-Nr.: 66147-1-Ig	GenBank-Zugangsnummer: BC107043	Reinigungsmethode: Protein-A-Reinigung
Größe: 150ul , Konzentration: 993 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): 27179	CloneNo.: 1F2B8
Wirt: Maus	Vollständiger Name: interleukin 1 family, member 6 (epsilon)	
Isotyp: IgG2b	Berechnete Masse: 158 aa, 18 kDa	
Immunogen Katalognummer: AG21487		

Anwendungen

Geprüfte Anwendungen:
WB,ELISA

Getestete Reaktivität:
Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Recombinant protein were subjected to SDS PAGE followed by western blot with 66147-1-Ig (IL-1F6 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.