

Allgemeine Informationen

Katalog-Nr.:

66312-1-Ig

Größe:

150ul, Konzentration: 1000 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Maus

Isotyp:

IgG2b

Immunogen Katalognummer:

AG24087

GenBank-Zugangsnummer:

NM-010552

GeneID (NCBI):

16171

Vollständiger Name:

interleukin 17A

Reinigungsmethode:

Protein-A-Reinigung

CloneNo.:

2F9D1

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

Il17a is a pro-inflammatory cytokine as a member of the interleukin-17 family. It plays a central role in host defense against diverse pathogens. It is produced by activated T-cells and certain cell types of innate immune system. The active protein functions as either a homodimer with other interleukin-17 family members and signals through the interleukin-17 receptor to induce inflammatory cytokine production. Aberrant expression of this gene is associated with autoinflammatory diseases including rheumatoid arthritis, psoriasis and multiple sclerosis.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

