

Nur für Forschungszwecke

THRA Monoklonaler Antikörper

Katalog-Nr.: 66703-1-Ig

Vorgestelltes Produkt



Allgemeine Informationen

Katalog-Nr.:

66703-1-Ig

Größe:

150ul

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG27189

GenBank-Zugangsnummer:

BC002728

GeneID (NCBI):

7067

Vollständiger Name:

thyroid hormone receptor, alpha
(erythroblastic leukemia viral (v-erb-
a) oncogene homolog, avian)

Berechnete Masse:

54 aa, 7 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

1A11C8

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



WB result of THRA antibody (66703-1-Ig; 1:4000; incubated at room temperature for 1.5 hours) with sh-Control and sh-THRA transfected A431 cells.