

Nur für Forschungszwecke

KIAA0284 Monoklonaler Antikörper

Katalog-Nr.:66749-1-Ig



Allgemeine Informationen

Katalog-Nr.:

66749-1-Ig

Größe:

150ul, Konzentration: 1800 µg/ml von283638

Nanodrop und 1000 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Maus

Isotyp:

IgG2b

Immunogen Katalognummer:

AG26852

GenBank-Zugangsnummer:

BC112928

GeneID (NCBI):

283638

Vollständiger Name:

KIAA0284

Berechnete Masse:

1589 aa, 172 kDa

Beobachtete Masse:

172 kDa

Reinigungsmethode:

Protein-A-Reinigung

CloneNo.:

1D4B8

Empfohlene Verdünnungen:

WB 1:1000-1:4000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : HeLa-Zellen, HEK-293-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 66749-1-Ig (KIAA0284 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.