

Allgemeine Informationen

Katalog-Nr.:

66962-1-Ig

Größe:

150ul, Konzentration: 1600 µg/ml von 5069

Nanodrop und 1000 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG11863

GenBank-Zugangsnummer:

BC078657

GeneID (NCBI):

1627

Vollständiger Name:

pregnancy-associated plasma protein A, pappalysin 1

Berechnete Masse:

1627 aa, 181 kDa

Beobachtete Masse:

180-220 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

1G5A7

Empfohlene Verdünnungen:

WB 1:4000-1:10000
IHC 1:500-1:1000

Anwendungen

Geprüfte Anwendungen:

IHC, WB, ELISA

In Publikationen genannte Anwendungen:

WB

Getestete Reaktivität:

Human

Zitierte Arten:

Human

Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.

Positivkontrollen:

WB : humanes Plazenta-Gewebe, MCF-7-Zellen

IHC : humanes Plazenta-Gewebe,

Hintergrundinformationen

Bemerkenswerte Veröffentlichungen

Verfasser	Pubmed ID	Journal	Anwendung
Go Itoh	34379869	Mol Oncol	WB

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




human placenta tissue were subjected to SDS PAGE followed by western blot with 66962-1-Ig (PAPPA antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.



WB result of PAPPA antibody (66962-1-Ig; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-PAPPA transfected MCF-7 cells.



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 66962-1-Ig (PAPPA antibody) at dilution of 1:1000 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 66962-1-Ig (PAPPA antibody) at dilution of 1:1000 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).