

Nur für Forschungszwecke

PRTN3 Monoklonaler Antikörper

Katalog-Nr.:67030-1-Ig



Allgemeine Informationen

Katalog-Nr.: 67030-1-Ig	GenBank-Zugangsnummer: BC096183	Reinigungsmethode: Protein-A-Reinigung
Größe: 150ul , Konzentration: 2100 µg/ml von5657 Nanodrop und 1000 µg/ml durch die Bradford-Methode mit BSA als Standard;	GeneID (NCBI): proteinase 3	CloneNo.: 2D1B8
Wirt: Maus	Vollständiger Name: proteinase 3	Empfohlene Verdünnungen: WB 1:1000-1:4000 IHC 1:1000-1:4000 IF 1:200-1:800
Isotyp: IgG2b	Berechnete Masse: 256 aa, 28 kDa	
Immunogen Katalognummer: AG25915	Beobachtete Masse: 28 kDa	

Anwendungen

Geprüfte Anwendungen: IF, IHC, WB, ELISA	Positivkontrollen:
Getestete Reaktivität: Human	WB : humanes Plazenta-Gewebe,
Hinweis-IHC: Antigendemaskierung mit TE-Puffer pH 9,0 empfohlen. (*) Wahlweise kann die Antigendemaskierung auch mit Citratpuffer pH 6,0 erfolgen.	IHC : humanes Milzgewebe, humanes Tonsillitisgewebe
	IF : humanes Tonsillitisgewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:
Bei -20°C lagern.
Lagerungspuffer:
PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Ausgewählte Validierungsdaten




human placenta tissue were subjected to SDS PAGE followed by western blot with 67030-1-Ig (PRTN3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human spleen tissue slide using 67030-1-Ig (PRTN3 antibody) at dilution of 1:2000 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human spleen tissue slide using 67030-1-Ig (PRTN3 antibody) at dilution of 1:2000 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using PRTN3 antibody (67030-1-Ig, Clone: 2D1B8) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).