

Nur für Forschungszwecke

DTX2 Monoklonaler Antikörper

Katalog-Nr.:67209-1-Ig



Allgemeine Informationen

Katalog-Nr.:

67209-1-Ig

Größe:

150ul, Konzentration: 2400 µg/ml von113878

Nanodrop und 1000 µg/ml durch die Bradford-Methode mit BSA als Standard;

Wirt:

Maus

Isotyp:

IgG2a

Immunogen Katalognummer:

AG12022

GenBank-Zugangsnummer:

BC064535

GeneID (NCBI):

113878

Vollständiger Name:

deltex homolog 2 (Drosophila)

Berechnete Masse:

253aa,28 kDa; 575aa,64 kDa

Beobachtete Masse:

68 kDa

Reinigungsmethode:

Protein-A-Reinigung

CloneNo.:

2D4F8

Empfohlene Verdünnungen:

WB 1:1000-1:6000

Anwendungen

Geprüfte Anwendungen:

WB,ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : HeLa-Zellen, A549-Zellen, HepG2-Zellen, MCF-7-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Various lysates were subjected to SDS PAGE followed by western blot with 67209-1-Ig (DTX2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.