

Nur für Forschungszwecke

MMP14 Monoklonaler Antikörper

Katalog-Nr.: 67225-1-Ig



Allgemeine Informationen

Katalog-Nr.:

67225-1-Ig

Größe:

150ul

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG6062

GenBank-Zugangsnummer:

BC064803

GeneID (NCBI):

4323

Vollständiger Name:

matrix metalloproteinase 14
(membrane-inserted)

Berechnete Masse:

66 kDa

Reinigungsmethode:

Protein-A-Reinigung

CloneNo.:

2F11A11

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

