

Nur für Forschungszwecke

ABTB1 Monoklonaler Antikörper

Katalog-Nr.: 67246-1-Ig



Allgemeine Informationen

Katalog-Nr.:

67246-1-Ig

Größe:

150ul, Konzentration: 1000 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG28959

GenBank-Zugangsnummer:

BC011858

GeneID (NCBI):

80325

Vollständiger Name:

ankyrin repeat and BTB (POZ) domain
containing 1

Berechnete Masse:

54 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

2D2B11

Empfohlene Verdünnungen:

WB 1:1000-1:4000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Hausschwein, Human

Positivkontrollen:

WB: Hausschwein-Cerebellum-Gewebe,

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



pig cerebellum tissue were subjected to SDS PAGE followed by western blot with 67246-1-Ig (ABTB1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.