

Nur für Forschungszwecke

BRCA2 Monoklonaler Antikörper

Katalog-Nr.: 67285-1-Ig

1 Publikationen



Allgemeine Informationen

Katalog-Nr.:

67285-1-Ig

Größe:

150ul, Konzentration: 553 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Maus

Isotyp:

IgG1

GenBank-Zugangsnummer:

NM_000059

GeneID (NCBI):

675

Vollständiger Name:

breast cancer 2, early onset

Berechnete Masse:

384 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

3G1C7

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

Bemerkenswerte Veröffentlichungen

Verfasser

Pubmed ID

Journal

Anwendung

Ishraq Alim

35616648

Genet Med

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

