

Nur für Forschungszwecke

ARL4C Monoklonaler Antikörper

Katalog-Nr.: 67356-1-Ig



Allgemeine Informationen

Katalog-Nr.:

67356-1-Ig

Größe:

150ul

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG29496

GenBank-Zugangsnummer:

BC001051

GeneID (NCBI):

10123

Vollständiger Name:

ADP-ribosylation factor-like 4C

Berechnete Masse:

21 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

1B8B5

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

