

Nur für Forschungszwecke

TRIM72 Monoklonaler Antikörper

Katalog-Nr.:67403-1-Ig



Allgemeine Informationen

Katalog-Nr.: 67403-1-Ig	GenBank-Zugangsnummer: BC033211	Reinigungsmethode: Protein-G-Reinigung
Größe: 150ul , Konzentration: 1300 µg/ml von493829 Nanodrop;	GeneID (NCBI): 493829	CloneNo.: 5F2C11
Wirt: Maus	Vollständiger Name: tripartite motif-containing 72	Empfohlene Verdünnungen: WB 1:1000-1:6000 IF 1:200-1:800
Isotyp: IgG1	Berechnete Masse: 477 aa, 53 kDa	
Immunogen Katalognummer: AG2560		

Anwendungen

Geprüfte Anwendungen: IF, WB, ELISA	Positivkontrollen:
Getestete Reaktivität: Hausschwein, Human, Maus	WB : Hausschwein-Skelettmuskelgewebe, Hausschwein-Herzgewebe
	IF : Maus-Skelettmuskelgewebe,

Hintergrundinformationen

Lagerung


Lagerungsbedingungen:
Bei -20°C lagern. Nach dem Versand ein Jahr lang stabil
Lagerungspuffer:
PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



pig heart tissue were subjected to SDS PAGE followed by western blot with 67403-1-Ig (TRIM72 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed mouse skeletal muscle tissue using TRIM72 antibody (67403-1-Ig, Clone: 5F2C11) at dilution of 1:400 and CoraLite@488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).



pig skeletal muscle tissue were subjected to SDS PAGE followed by western blot with 67403-1-Ig (TRIM72 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed mouse skeletal muscle tissue using TRIM72 antibody (67403-1-Ig, Clone: 5F2C11) at dilution of 1:400 and CoraLite@488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).