

SET Monoklonaler Antikörper

Katalog-Nr.: 67478-1-Ig

Allgemeine Informationen

Katalog-Nr.:

67478-1-Ig

GenBank-Zugangsnummer:

NM_003011

Reinigungsmethode:

Protein-A-Reinigung

Größe:

150ul, Konzentration: 2200 µg/ml von 6418

GeneID (NCBI):

6418

CloneNo.:

2B4D1

Nanodrop und 1000 µg/ml durch die Bradford-Methode mit BSA als Standard;

Vollständiger Name:

SET nuclear oncogene

Empfohlene Verdünnungen:

WB 1:1000-1:6000

Wirt:

Maus

Berechnete Masse:

33 kDa

Isotyp:

IgG2b

Beobachtete Masse:

33-36 kDa

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Positivkontrollen:

WB: Jurkat-Zellen, 4T1-Zellen, HeLa-Zellen, HepG2-Zellen

Getestete Reaktivität:

Human, Maus, Ratte

Hintergrundinformationen

SET, also named as IGAAD, I-2PP2A, TAF-I, IGAAD and PHAPII, belongs to the nucleosome assembly protein (NAP) family. SET is a multifunctional oncoprotein which is ubiquitously expressed in all kinds of cells. The SET protein participates in many cellular processes including cell cycle, cell migration, apoptosis, transcription, and DNA repair. The SET protein plays important roles in facilitating cell proliferation via many pathways that promote tumorigenesis and metastasis.

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS mit 0.02% Natriumazid und 50% Glycerin pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Ausgewählte Validierungsdaten



Jurkat cells were subjected to SDS PAGE followed by western blot with 67478-1-Ig (SET antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.