

Allgemeine Informationen

Katalog-Nr.:

67718-1-Ig

Größe:

150ul, Konzentration: 966 µg/ml
durch die Bradford-Methode mit BSA
als Standard;

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG30382

GenBank-Zugangsnummer:

BC014949

GeneID (NCBI):

79132

Vollständiger Name:

DEXH (Asp-Glu-X-His) box
polypeptide 58

Berechnete Masse:

678 aa, 77 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

2H1A1

Anwendungen

Geprüfte Anwendungen:

ELISA

Getestete Reaktivität:

Human

Hintergrundinformationen

DHX58 (Probable ATP-dependent RNA helicase DHX58) is also named as D11LGP2E, LGP2 and belongs to the RLR subfamily. DHX58, originally identified as a highly expressed gene in mammary tumors, is another cytoplasmic DEX(D/H)-box helicase that can recognize RNA (PMID: 18411269). It acts as a positive, but not negative, regulator of RIG-I- and MDA5-dependent recognition of RNA virus infection and plays a pivotal role in antiviral responses in vivo (PMID: 20080593).

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

***** 20ul-Größen enthalten 0.1% BSA**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

