

Nur für Forschungszwecke

GIMAP2 Monoklonaler Antikörper

Katalog-Nr.: 67720-1-Ig



Allgemeine Informationen

Katalog-Nr.:
67720-1-Ig

Größe:

150ul, Konzentration: 1500 µg/ml von 26157

Nanodrop und 749 µg/ml durch die
Bradford-Methode mit BSA als
Standard;

Wirt:

Maus

Isotyp:

IgG1

Immunogen Katalognummer:

AG4368

GenBank-Zugangsnummer:

BC032345

GeneID (NCBI):

von 26157

Vollständiger Name:

GTPase, IMAF family member 2

Berechnete Masse:

337 aa, 38 kDa

Beobachtete Masse:

33-38 kDa

Reinigungsmethode:

Protein-G-Reinigung

CloneNo.:

1F2A2

Empfohlene Verdünnungen:

WB 1:1000-1:4000

Anwendungen

Geprüfte Anwendungen:

WB, ELISA

Getestete Reaktivität:

Human

Positivkontrollen:

WB : MOLT-4-Zellen, Jurkat-Zellen

Hintergrundinformationen

Lagerung

Lagerungsbedingungen:

Bei -20°C lagern.

Lagerungspuffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Aliquotieren ist nicht notwendig bei -20°C Lagerung

*** 20ul-Größen enthalten 0.1% BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Ausgewählte Validierungsdaten



MOLT-4 cells were subjected to SDS PAGE followed by western blot with 67720-1-Ig (GIMAP2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.